

说明书尺寸:110\*80mm

80 mm

110 mm

habo

OGAWA

User Manual

Please read this manual carefully before use.

Wireless Charging Blackhead Remover

Dear customers, thank you for choosing wireless charging blackhead remover. Please read the instruction carefully and operate correctly in accordance with the operating manual before using the product.

Product Structure:

Probes

Blackheads Removing/Lighting On/Off Button

Gears

Gears Adjusting Button

LED Light

Indicator Light of Charging

Wireless Charging Induction Area

USB Charging Port

Microcrystalline Probe

Oval Probe

Small Round Probe

Big Round Probe

Product Specifications:

Rated Current: 200mA	Charging Voltage: 5V/1A
Power Supply: Lithium battery	Rated Output Power: 0.75W
Battery Voltage: 3.7V	
Charging Method: Wireless charging	

Features:

Purify Pores: Deeply remove blackheads, acnes and facial oil inside skin and pores.

Product Instruction

Using Methods

1. Clean your face, and you can also use the blackhead remover and facial steamer to open the pore. The closed pores will reduce the cleaning effect.

closed pores will reduce the cleaning effect.  
For the first time use, it is recommended to test the suction with your forearm before apply it to your face.  
In order to avoid the injury, do not let the probe stay in one place for too long, and do not press it hardly or pull vertically. The improper use may cause facial redness, which will fade away in few days.  
Please do not apply to the skin wounds or scars.  
Please keep the body away from water or other liquid.  
Please do not use chemical cleaner so as to avoid corrosion.  
It is recommended that the cleaning should be once a week for dry skin, twice for oily skin, and each time should be within 5 minutes. If it is too frequently, it may lead to skin allergies.  
Please clean the probes every time after using to avoid the filter blocked. If you plan to use the alcohol, the concentration should be less than 70%.  
Please do not use the product on children in case of damaging the skin.  
This product is equipped with a USB charging cable, but you may need to purchase an additional 5V/1A adapter to work with it.

Package Contents

1 x Wireless Charging Blackhead Remover  
1 x Big Round Probe  
1 x Small Round Probe  
1 x Oval Probe  
1 x Microcrystalline Probe  
1 x Wireless Charging Base  
1 x Charging Cable  
1 x Pack of Isolation Sponge Pads  
1 x Stainless Steel Blackhead Removal Needle  
1 x User Manual

ContactUs

MALAYSIA  
Healthy World Lifestyle Sdn. Bhd  
No.22 Jalan Anggerik Mokara 31/47,  
Kota Kemuning, 40460 Shah Alam,  
Selangor Darul Ehsan, Malaysia  
Tel: +603-5121 4286  
Fax: +603-5121 4386

2. Insert the isolation pads into the slot at the top of the cleaner, select the proper probe and install it on the nozzle.

3. Short press the ON/OFF button to open the cleaner. The default is the Soft level (Our product will automatically memorize your selected level before shutdown in the following uses).

4. According to personal need, short press up or down of the Levels Adjusting Button to choose the proper level. For the first time use, it is recommended to test the suction with your forearm before apply it to your face.

5. Gently hold the probe close to your face to remove blackheads, and slowly move the probe. In order to avoid the injury, do not let the probe stay in one place for too long, and do not press it hardly or pull vertically.

6. After cleaning, short press the ON/OFF button to shut down the system.  
7. Remove the probe, wipe the dirt off with cotton swab, then wash it clean, and you can dry it with a towel or let it dry naturally. Please avoid to wash or immerse the body in water.

8. If the natural light outside is insufficient, you are able to long press ON/OFF button to turn on or shutdown LED lighting system.

Charging Method

1. When the battery is low, the red indicator flashes and goes out after several seconds to remind you to charge the device.  
2. Connect the charging cable with a 5V/1A adapter (not included) and the wireless charging base.  
3. Put the blackhead remover on the wireless charging base, the electric quantity indicator keeps red light on.  
4. When fully charged, the electric quantity indicator keeps blue light on instead, then just plug the charging cable out (Charging time is about 8 hours).

Notes

Due to the long distance transportation, the battery may be out of power. If it does not work after starting up, it is recommended to use the product after a full charge.  
To ensure a better cleaning, please use the blackhead remover or facial steamer before use the product to open the pore. The